

## ABSTRACT

In support of NASA ARMD's code validation project, we have made significant progress by providing the first quantitative single-shot multi-scalar data from a turbulent elevated-pressure (5 atm), swirl-stabilized, lean direct injection (LDI) type research burner operating on CH<sub>4</sub>-air using a spatially-resolved pulsed-laser spontaneous Raman diagnostic technique. The Raman diagnostics apparatus and data analysis that we present here were developed over the past 6 years at Glenn Research Center. From the Raman scattering data, we produce spatially-mapped probability density functions (PDF's) of the instantaneous temperature, determined using a newly developed low-resolution effective rotational bandwidth (ERB) technique. The measured 3-scalar (triplet) correlations, between temperature, CH<sub>4</sub>, and O<sub>2</sub> concentrations, as well as their PDF's, also provide a high-level of detail into the nature and extent of the turbulent mixing process and its impact on chemical reactions in a realistic gas turbine injector flame at elevated pressures. The multi-scalar triplet data set presented here provides a good validation case for CFD combustion codes to simulate by providing both average and statistical values for the 3 measured scalars.



# **Single-Shot Scalar-Triplet Measurements in High-Pressure Swirl-Stabilized Flames for Combustion Code Validation**

**Jun Kojima\***  
*Ohio Aerospace Institute*

**Quang-Viet Nguyen**  
*NASA GRC (RTB)*

HPGB Facility ERB/SE-5  
H<sub>2</sub>-Air Flame at 20 atm


\* Principal Investigator  
FY06 NRA


**Goal:** Provide critical chemical and physical data in turbulent combustion at realistic subsonic/supersonic cruise condition for **validating predictive low-emissions combustor codes**

**Challenges:** Successful accurate multiscale measurements (temperature, major species mole fractions, mixture fraction, velocity) in optically-harsh droplet-laden flows at **high pressures**.

### Strategy (tasks):


- ✓ Develop a quantitative time-resolved **laser Raman diagnostics**
- ✓ A series of experiments of increasing flow complexity: (i) Calibration burner (e.g., flat-flame burner), (ii) Realistic concept burners (swirl-stabilized burner) at **elevated pressures** up to 10~30 atm
- ✓ Address the effects of chemical complexity:  $H_2$ ,  $CH_4$ , single-component liquid (e.g., hexane, iso-octane, n-heptane), jet-A fuel





- Temperature and species reference
- Diagnostics calibration
- Gaseous fuels
- Pre-vaporized single-component liquid fuel (electronically controlled fuel vaporizer and heating system)
- Support visible and UV laser Raman systems
- Study on optical characteristics of sooting flames and liquid fuel combustion chemistry
- Support the Aeronautics milestones to mitigate risks due to PSO holding up high-pressure experiments



Flat-flame burner



Fuel/air mixture  
(gas, or pre-vaporized liquid fuel)





## High Pressure Burner Facility (SE-5)



- ✓ Remotely controlled with auto process controller
- ✓ Pressure up to **30 atm** (currently 10 atm requested for safety permit)
- ✓ **Versatile** for burner platform (calibration burner, turbulent jet, LDI)
- ✓ Optical access (4 ports)
- ✓ Gaseous and **liquid fuel** capacity
- ✓ Air **pre-heater** installed (up to 1200F)

# Raman Scattering Diagnostic System




## Raman Spectra Observed in Combustion

### Advantage

- Raman scattering –simultaneous multiple species concentration and temperature
- Time-average (mean) or Single-shot (instantaneous; rms)
- Quantitative


### Challenge

- Spectral interference ('Cross-talk')




# Raman Calibration Experiments & Simulation


## Platform



Static cell



Calib. burner



Raman spectrum  
simulation

## Fuel (lean ~ rich)

H <sub>2</sub>
H <sub>2</sub> -CO
CH <sub>4</sub>
Single-component liquid (pre-vaporized)
Single-component liquid
Multi-component (jet-A)

Air
-----

## Raman Instrumentation

Pump lasers: 532/355/266 nm
Thermometry: Vib. / Rot.
Polarization property




## Matrix elements, $k_{i,j}(T)$ (Cross-talks)

Major species (H <sub>2</sub> , O <sub>2</sub> , N <sub>2</sub> , CO, CO <sub>2</sub> , H <sub>2</sub> O, HC's)
H <sub>2</sub> ► CO
H <sub>2</sub> ► CO <sub>2</sub>
H <sub>2</sub> ► H <sub>2</sub> O
N <sub>2</sub> ► CO
O <sub>2</sub> ► CO <sub>2</sub>
CO <sub>2</sub> ► N <sub>2</sub>
C <sub>2</sub> * ► N <sub>2</sub> , CO, HC's
PAH ► Broadband background



# Determine the 'Cross-talk' Calibration Matrix

## Raman Spectra Library





Calibration  
matrix

— High temp.  
— Low temp.

Mole frac.


$$N_i = (1/E_{\text{laser}}) k_{i,j}(T) S_i \quad \text{Raman signal}$$

## Calibrated Results




# NASA Lean Direct Injection (LDI) Swirl-Stabilized Research Burner — Preliminary Test —

Photos of LDI flame and burner face




Single Gaseous LDI Injector Design  
(Solid model for CFD)



- Swirl-stabilized direct fuel injection design (gaseous)
- Integrated with existing high-pressure rig
- 6 jets (0.8 mm in dia.) angled at 45 deg to burner axis
- Initial test on H<sub>2</sub> and CH<sub>4</sub> fuels with unheated air at 5 atm
- Collaborating with National Combustor Code (NCC)

## Single-Shot (Time-Resolved) Raman Measurements


Time variation of single-shot Raman spectra  
in CH<sub>4</sub>-air LDI flame  
(5 atm,  $\Phi = 0.5$ )






- Temperature from the rotational spectra has been measured with accuracy 7%
- Data permits **statistical PDF's** of temperature and species
- Provides a signature for **characterizing degree of mixing and reaction**

- One shot = One instance (space-time point)
- Single-shot Raman data shows “random” change due to turbulence
- Direct and simultaneous measurement of fuel/oxidizer concentrations and temperature (with data processing with calibration matrix)
- Developed **new thermometry approach** (rot. bandwidth) with high SNR

Rotational N<sub>2</sub> Raman band




## Processing the Single-Shot Data



## Probability Density Functions: PDF's (Temperature, CH<sub>4</sub>, and O<sub>2</sub>)

CH<sub>4</sub>-air LDI flame (5 atm,  $\Phi = 0.5$ )





✓ Spatially-mapped multi-scalar PDF “library” in given conditions provides excellent *validation data* for a predictive combustion code.

✓ Difference between mean and distribution (statistical value)



# Multiscalar Analysis in 5-atm CH<sub>4</sub>-Air LDI Flame





*'Production-mode' data*  
: Direct output from a MATLAB code

1. Temperature PDF's
2. Temp. vs species correlations  
(red: O<sub>2</sub>; blue: CH<sub>4</sub>)
3. Triplet correlations (temp-ox-fuel)  
(with a global Phi of 0.5 line)

## Data interpretation

- ✓ Hot spots (high NO<sub>x</sub>)
- ✓ Unique bi-modal PDF's — recirculation zone, or combustion oscillations
- ✓ Turbulent-chemistry interaction
- ✓ Fuel-air mixing characteristics
- ✓ Unburnt pockets

# Multiscalar Analysis in 5-atm CH<sub>4</sub>-Air LDI Flame



## Data interpretation

- ✓ Gaussian-like narrow Temp. distribution (centered around the adiabatic temp. at Phi of 0.5)
- ✓ No fuel (CH<sub>4</sub>) residual = fully consumed
- ✓ No hot or cold spots



✓ Homogeneous, well-reacted post-flame zone

## Conclusions



- A single-shot (time-resolve) capability of a laser Raman diagnostics has been confirmed.
- Preliminary data of time-resolved multispectral data in a high-pressure (turbulent) swirl-stabilized flame has been acquired.
- A new single-shot data-processing scheme (computer code) has been developed for 'production mode' thermo-chemical analysis.
- Scalar PDF's and 3D (temp-oxy-fuel) correlations showed promising capability of future use in code validation

## Work-in-progress

- APCD (atmospheric pressure combustion diagnostics) facility is under construction to calibrate and improve the Raman diagnostics applicable to liquid fuels
- Modified visible and UV Raman diagnostic systems is under development to cope with harsh environments
- Computer code to simulate Raman spectra of major species including CO<sub>2</sub> and H<sub>2</sub>O (except hydrocarbon) is under development to complete the calibration matrix
- Pressure Safety Office safety permits, variances in process.